



# Polyurethane Peripherally Inserted Central Catheter with Peel-Away Introducer Sheaths



The hub is clearly marked with "PICC" and French size to readily identify catheter type.

The Footprint® introducer sheaths have unique design features to enhance first stick success.



## Peripherally Inserted Central Catheter

- Designed especially for neonatal/pediatric patients
- Polyurethane catheter with integrated extension set
- Soft stabilization wing for enhanced securement
- Indelible markings every centimeter
- Fully radiopaque
- Latex and DEHP free

## Introducer Sheath

- Notched needle with minimal lie distance for immediate visualization of flashback designed to increase first stick success
- Low profile, compact introducer designed to enhance small vein access



**Footprint**<sup>®</sup>  
medical inc.

210.226.2600

[www.footprintmed.com](http://www.footprintmed.com)



# Ordering Information:

## Polyurethane Peripherally Inserted Central Catheters

P1PIC1.4-S	Polyurethane PICC Complete Kit <i>(w/ single-lumen catheter &amp; introducer)</i>	1.4 Fr.	5/cs
P1PIC1.4-C	Polyurethane PICC Single-Lumen Catheter Only	1.4 Fr.	5/cs
P1PIC1.9-S	Polyurethane PICC Complete Kit <i>(w/ single-lumen catheter &amp; introducer)</i>	1.9 Fr.	5/cs
P1PIC1.9-C	Polyurethane PICC Single-Lumen Catheter Only	1.9 Fr.	5/cs
P2PIC1.9-S	Polyurethane PICC Complete Kit <i>(w/ dual-lumen catheter &amp; introducer)</i>	1.9 Fr.	5/cs
P2PIC1.9-C	Polyurethane PICC Dual-Lumen Catheter Only	1.9 Fr.	5/cs
DC-PIC	PICC Dressing Change Kit		10/cs

## Peel-Away Introducer Sheaths

SHEATH1.4	PICC Sheath Introducer Only	1.4 Fr.	5/cs
SHEATH1.9	PICC Sheath Introducer Only	1.9 Fr.	5/cs

*"Helping tiny babies make their own footprints in the world..."*



12727 Cimarron Path  
San Antonio, TX 78249  
210.226.2600  
[www.footprintmed.com](http://www.footprintmed.com)