



# Silicone Peripherally Inserted Central Catheter with Peel-Away Introducer Sheaths



The Footprint® introducer sheaths have unique design features to enhance first stick success.



The hub is clearly marked with "PICC" and French size to readily identify catheter type.

## Peripherally Inserted Central Catheter

- Designed especially for neonatal/pediatric patients
- Silicone catheter with integrated extension set
- Soft stabilization wing for enhanced securement
- Indelible markings every centimeter
- Fully radiopaque
- Latex and DEHP free

## Introducer Sheath

- Notched needle with minimal lie distance for immediate visualization of flashback designed to increase first stick success
- Low profile, compact introducer designed to enhance small vein access



**Footprint**<sup>®</sup>  
medical inc.

210.226.2600  
www.footprintmed.com



## Ordering Information:

### Silicone Peripherally Inserted Central Catheters

S1PIC1.9-S	Silicone PICC Complete Kit <i>(with catheter &amp; introducer)</i>	1.9 Fr.	5/cs
S1PIC1.9-C	Silicone PICC Catheter Only	1.9 Fr.	5/cs
S1PIC1.9-L	Silicone PICC Catheter Only (43cm)	1.9 Fr.	5/cs
DC-PIC	PICC Dressing Change Kit		10/cs

### Peel-Away Introducer Sheaths

SHEATH1.9	PICC Sheath Introducer Only	1.9 Fr.	5/cs
-----------	-----------------------------	---------	------

*"Helping tiny babies make their own footprints in the world..."*



12727 Cimarron Path  
San Antonio, TX 78249  
210.226.2600  
[www.footprintmed.com](http://www.footprintmed.com)