



# Urinary Drainage Catheter



The hub is clearly marked with "URINE" and French size to readily identify the catheter type.



The 100ml burette provides accurate measurement of urinary output.

## Urinary Drainage Catheter

- Made of soft, biocompatible medical grade silicone designed to enhance extended dwell protocols
- Bold, indelible markings every centimeter
- Catheter has a rounded tip with two precision drilled side holes to reduce trauma during insertion and minimize urethral irritation
- Luer-Lok™ hub with integrated closure plug
- Latex and DEHP free with radiopaque stripe
- Also available in Urine Collection Kit
- Closed System Kit Contains:
  - Silicone urinary catheter
  - 100ml Burette
  - Antiseptic swab sticks
  - Sterile catheter lubricant
  - Fenestrated drape
  - Poly-lined drape
  - Gauze pads



**Footprint**  
medical inc.

210.226.2600  
www.footprintmed.com

Ordering Information for:

# Urinary Drainage Catheter



## Silicone Urinary Drainage Catheters

S1UDC3.5	3.5 Silicone Urinary Drainage Catheter	10/cs
S1UDC5.0	5.0 Silicone Urinary Drainage Catheter	10/cs
S1UDC6.5	6.5 Silicone Urinary Drainage Catheter	10/cs
S1UDC8.0	8.0 Silicone Urinary Drainage Catheter	10/cs

## Silicone Urinary Drainage Kit without Catheter

UDKNC	Closed System Urinary Collection Kit without Catheter	5/cs
-------	---	------

## Silicone Urinary Collection Kit with Catheter

S1UDK3.5	3.5 Fr. Silicone Urinary Collection Kit with Catheter	5/cs
S1UDK5.0	5.0 Fr. Silicone Urinary Collection Kit with Catheter	5/cs
S1UDK6.5	6.5 Fr. Silicone Urinary Collection Kit with Catheter	5/cs
S1UDK8.0	8.0 Fr. Silicone Urinary Collection Kit with Catheter	5/cs

*"Helping tiny babies make their own footprints in the world..."*



12727 Cimarron Path  
San Antonio, TX 78249  
210.226.2600  
[www.footprintmed.com](http://www.footprintmed.com)